

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TEKT4 Polyclonal Antibody

多克隆抗体

产品货号: PA07878

第三版

**描述:**Structural component of ciliary and flagellar microtubules. Forms filamentous polymers in the walls of ciliary and flagellar microtubules. Found in the abaxial (convex) surface of outer dense fibers in sperm flagella. TEKT4 is necessary for the proper coordinated beating of the sperm flagellum and male reproductive physiology.

**抗原:**Recombinant protein of human TEKT4

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**