

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TRPC7 Polyclonal Antibody

多克隆抗体

产品货号: PA07930

第三版

**描述:** This protein may form a receptor-activated non-selective calcium permeant cation channel. Probably is operated by a phosphatidylinositol second messenger system activated by receptor tyrosine kinases or G-protein coupled receptors. Activated by diacylglycerol (DAG) . May also be activated by intracellular calcium store depletion. It can interact with Human, Mouse X1 and RNF24. Interacts (via ANK-repeat domains) with PRKG1.

**抗原:** Recombinant protein of human TRPC7

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**