

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TSHR Polyclonal Antibody

多克隆抗体

产品货号: PA07945

第三版

**描述:**The protein encoded by this gene is a membrane protein and a major controller of thyroid cell metabolism. The encoded protein is a receptor for thyrothopin and thyrostimulin, and its activity is mediated by adenylate cyclase. Defects in this gene are a cause of several types of hyperthyroidism. Three transcript variants encoding different isoforms have been found for this gene.

**抗原:**Recombinant protein of human TSHR

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$  g/ml.

**特异性:**Expressed in the thyroid.