

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BPIFB1 Polyclonal Antibody

多克隆抗体

产品货号: PA07961

第三版

描述:The protein encoded by this gene may be involved in the innate immune response to bacterial exposure in the mouth, nasal cavities, and lungs. The encoded protein is secreted and is a member of the BPI/LBP/PLUNC protein superfamily. This gene is found with other members of the superfamily in a cluster on chromosome 20.

抗原:Synthetic peptide of human BPIFB1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in trachea, nasal septal epithelium and lung.