

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SERPINA5 Polyclonal Antibody

多克隆抗体

产品货号: PA07974

第三版

**描述:**The protein encoded by this gene is a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. This gene is one in a cluster of serpin genes located on the q arm of chromosome 14. This family member is a glycoprotein that can inhibit several serine proteases, including protein C and various plasminogen activators and kallikreins, and it thus plays diverse roles in hemostasis and thrombosis in multiple organs.

**抗原:**Recombinant protein of human SERPINA5

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed by the liver and secreted in plasma.