

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human STK36 Polyclonal Antibody

多克隆抗体

产品货号: PA07975

第三版

描述: This gene encodes a member of the serine/threonine kinase family of enzymes. This family member is similar to a *Drosophila* protein that plays a key role in the Hedgehog signaling pathway. This human protein is a positive regulator of the GLI zinc-finger transcription factors. Knockout studies of the homologous mouse gene suggest that defects in this human gene may lead to congenital hydrocephalus, possibly due to a functional defect in motile cilia.

抗原: Synthetic peptide of human STK36

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed at low levels in most fetal tissues, adult ovaries and at high levels in adult testis, where it is localized in germ cells.1
Publication