

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SMC6 Polyclonal Antibody

多克隆抗体

产品货号: PA08054

第三版

描述:Structural maintenance of chromosomes protein 6 is a protein that in humans is encoded by the SMC6 gene. Core component of the SMC5-SMC6 complex, a complex involved in DNA double-strand breaks by homologous recombination. The complex may promote sister chromatid homologous recombination by recruiting the SMC1-SMC3 cohesin complex to double-strand breaks.

抗原:Synthetic peptide of human SMC6

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed. Overexpressed in testis.