

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ACOT1 Polyclonal Antibody

多克隆抗体

产品货号: PA08057

第三版

描述: Acyl-CoA thioesterases such as ACOT1, hydrolyze acyl-CoAs to the free fatty acid and CoA. ACOTs therefore play key roles in maintaining the intracellular ratio between CoA esters of various lipids and free fatty acids. Acyl-CoA thioesterases are a group of enzymes that catalyze the hydrolysis of acyl-CoAs to the free fatty acid and coenzyme A (CoASH), providing the potential to regulate intracellular levels of acyl-CoAs, free fatty acids and CoASH. Active towards fatty acyl-CoA with chain-lengths of C12-C16.

抗原: Recombinant protein of human ACOT1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: