

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ACSM5 Polyclonal Antibody

多克隆抗体

产品货号: PA08075

第三版

描述:Fatty acids are incorporated into membranes and signaling molecules and have roles in energy storage and metabolism. These essential functions require activation of the fatty acid by acyl-coenzyme A (CoA) synthetases, such as ACSM5, which form an activating thioester linkage between the fatty acid and CoA. Has medium-chain fatty acid:CoA ligase activity with broad substrate specificity (in vitro). Acts on acids from C4 to C11 and on the corresponding 3-hydroxy- and 2,3- or 3,4-unsaturated acids (in vitro).

抗原:Recombinant protein of human ACSM5

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in kidney and liver.