

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPR38 Polyclonal Antibody

多克隆抗体

产品货号: PA08079

第三版

**描述:** Motilin is a 22 amino acid peptide hormone expressed throughout the gastrointestinal (GI) tract. The protein encoded by this gene is a motilin receptor which is a member of the G-protein coupled receptor 1 family. This member is a multi-pass transmembrane protein, and is an important therapeutic target for the treatment of hypomotility disorders.

**抗原:** Synthetic peptide of human MLNR

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expressed only in thyroid, stomach, and bone marrow.