

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ACTL8 Polyclonal Antibody

多克隆抗体

产品货号: PA08084

第三版

**描述:** This protein belongs to the actin family. Actin-like 8 is a protein in humans that is encoded by the ACTL8 gene. The high levels expression in testis and pancreas and weak expression in placenta. The subcellular location of this protein in cytoplasm, cytoskeleton.

**抗原:** Recombinant protein of human ACTL8

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Strongly expressed in testis and pancreas. Weak expression in placenta.