

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ACY3 Polyclonal Antibody

多克隆抗体

产品货号: PA08096

第三版

**描述:** Acy3, which encodes a deduced 318-amino acid protein with a calculated molecular mass of about 35.2 kD. Acy3 has no membrane domain, but it contains consensus sites for N-glycosylation, tyrosine sulfation, tyrosine and threonine phosphorylation, and myristoylation. Northern blot analysis detected high expression of an approximately 1.4-kb Acy3 transcript in kidney, followed by liver. Weaker expression was detected in heart, small intestine, brain, lung, testis, and stomach. Testis expressed an additional transcript of about 2.2 kb.

**抗原:** Recombinant protein of human ACY3

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:**