

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SLC43A2 Polyclonal Antibody

多克隆抗体

产品货号: PA08110

第三版

描述:System L amino acid transporters, such as SLC43A2, mediate sodium-independent transport of bulky neutral amino acids across cell membranes.

抗原:Synthetic peptide of human SLC43A2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in several tissues with, higher expression in placenta, kidney and peripheral blood leukocytes. In the kidney is detected in epithelial cells of the distal tubule and collecting duct. In the intestine is expressed mainly in crypt cells of the intestinal microvilli and epithelial cells in the base of the villus.