

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ADAM7 Polyclonal Antibody

多克隆抗体

产品货号: PA08116

第三版

描述: This gene encodes a member of the ADAMs family of zinc proteases. These transmembrane proteins play roles in multiple processes including cell signaling, adhesion and migration. The encoded protein lacks protease activity and may play roles in protein-protein interactions and cell adhesion processes including sperm-egg fusion. Mutations in this gene may be involved in the progression of melanoma.

抗原: Recombinant protein of human ADAM7

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Not detected in healthy melanocytes. Expressed in melanoma cells.