
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC16A12 Polyclonal Antibody

多克隆抗体

产品货号: PA08151

第三版

描述: This gene encodes a transmembrane transporter that likely plays a role in monocarboxylic acid transport. A mutation in this gene has been associated with juvenile cataracts with microcornea and renal glucosuria.

抗原: Synthetic peptide of human SLC16A12

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Most highly expressed in kidney, followed by retina, lung, and testis. Very weakly expressed in brain and liver. Also detected in lens.