

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC16A9 Polyclonal Antibody

多克隆抗体

产品货号: PA08157

第三版

描述: MCT9 (monocarboxylate transporter 9), also known as SLC16A9 (solute carrier family 16, member 9 (monocarboxylic acid transporter 9)) or C10orf36, is a 509 amino acid multi-pass membrane protein that belongs to the major facilitator superfamily. Like other MCT proteins, MCT9 functions as a monocarboxylate transporter that catalyzes the rapid proton-linked transport of monocarboxylates across the plasma membrane.

抗原: Synthetic peptide of human SLC16A9

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: