

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human P2RY2 Polyclonal Antibody

多克隆抗体

产品货号: PA08159

第三版

描述:The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor is responsive to both adenosine and uridine nucleotides. It may participate in control of the cell cycle of endometrial carcinoma cells. Three transcript variants encoding the same protein have been identified for this gene.

抗原:Synthetic peptide of human P2RY2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Spleen, testis, kidney, liver, lung, heart and brain.