

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SIGLEC11 Polyclonal Antibody

多克隆抗体

产品货号: PA08176

第三版

**描述:** This gene encodes a member of the sialic acid-binding immunoglobulin-like lectin family. These cell surface lectins are characterized by structural motifs in the immunoglobulin (Ig)-like domains and sialic acid recognition sites in the first Ig V set domain. This family member mediates anti-inflammatory and immunosuppressive signaling. Multiple transcript variants encoding different isoforms have been found for this gene.

**抗原:** Synthetic peptide of human SIGLEC11

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Expressed by macrophages in various tissues including Kupffer cells. Also found in brain microglia.