

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NES Polyclonal Antibody

多克隆抗体

产品货号: PA08182

第三版

**描述:** This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells. Required for brain and eye development. Promotes the disassembly of phosphorylated vimentin intermediate filaments (IF) during mitosis and may play a role in the trafficking and distribution of IF proteins and other cellular factors to daughter cells during progenitor cell division.

**抗原:** Synthetic peptide of human NES

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** CNS stem cells.