

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GFRA4 Polyclonal Antibody

多克隆抗体

产品货号: PA08300

第三版

描述:The protein encoded by this gene is a member of the GDNF receptor family. It is a glycosylphosphatidylinositol(GPI)-linked cell surface receptor for persephin, and mediates activation of the RET tyrosine kinase receptor. This gene is a candidate gene for RET-associated diseases. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

抗原:Synthetic peptide of human GFRA4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Predominantly expressed in the adult thyroid gland. Low levels also found in fetal adrenal and thyroid glands.