
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ANKRD54 Polyclonal Antibody

多克隆抗体

产品货号: PA08306

第三版

描述: Plays an important role in regulating intracellular signaling events associated with erythroid terminal differentiation. Interacts (via ankyrin repeat region) with LYN (via SH3-domain) in an activation-independent status of LYN By similarity. Forms a multiprotein complex with LYN and HCLS1 By similarity. Interacts with TSN2, VAV1, DBNL AND LASP1.

抗原: Recombinant protein of human ANKRD54

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: