

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat TPM2 Polyclonal Antibody

多克隆抗体

产品货号: PA08331

第三版

**描述:**This gene encodes beta-tropomyosin, a member of the actin filament binding protein family, and mainly expressed in slow, type 1 muscle fibers. Mutations in this gene can alter the expression of other sarcomeric tropomyosin proteins, and cause cap disease, nemaline myopathy and distal arthrogyrosis syndromes. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**抗原:**Synthetic peptide of human TPM2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Present in primary breast cancer tissue, absent from normal breast tissue.