

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NOS2 Polyclonal Antibody

多克隆抗体

产品货号: PA08341

第三版

描述: Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. This gene encodes a nitric oxide synthase which is expressed in liver and is inducible by a combination of lipopolysaccharide and certain cytokines. Three related pseudogenes are located within the Smith-Magenis syndrome region on chromosome 17.

抗原: Synthetic peptide of human NOS2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in the liver, retina, bone cells and airway epithelial cells of the lung. Not expressed in the platelets.