

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ATG12 Polyclonal Antibody

多克隆抗体

产品货号: PA08342

第三版

描述:Autophagy is a process of bulk protein degradation in which cytoplasmic components, including organelles, are enclosed in double-membrane structures called autophagosomes and delivered to lysosomes or vacuoles for degradation. ATG12 is the human homolog of a yeast protein involved in autophagy.

抗原:Recombinant protein of human ATG12

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.