
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat APTX Polyclonal Antibody

多克隆抗体

产品货号: PA08352

第三版

描述:This gene encodes a member of the histidine triad (HIT) superfamily. The encoded protein may play a role in single-stranded DNA repair through its nucleotide-binding activity and its diadenosine polyphosphate hydrolase activity. Mutations in this gene have been associated with ataxia-ocular apraxia. Alternatively spliced transcript variants have been identified for this gene.

抗原:Recombinant protein of human APTX

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed. In brain, it is expressed in the posterior cortex, cerebellum, hippocampus and olfactory bulb. Isoform 1 is highly expressed in the cerebral cortex and cerebellum, compared to isoform 2.