

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse DOPEY1 Polyclonal Antibody

多克隆抗体

产品货号: PA08356

第三版

**描述:** This gene belongs to the dopey family. It is a transporter protein, possibly involved in protein traffic between late Golgi and early endosomes. This gene plays a potential role in functional brain alterations and in the pathogenesis of mental retardation in Down syndrome, its overexpression in the brain regions, that are altered in Down syndrome patients and involved in learning and memory processes.

**抗原:** Synthetic peptide of human DOPEY1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4)，1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**