

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CD71 Polyclonal Antibody

多克隆抗体

产品货号: PA08368

第三版

描述:Transferrin receptor protein 1 (TfR1) also known as (Cluster of Differentiation 71) (CD71) is a protein that in humans is encoded by the TFRC gene. TfR1 is required for iron delivery from transferrin to cells. TfR1 is a transmembrane glycoprotein composed of two disulfide-linked monomers joined by two disulfide bonds. Each monomer binds one holo-transferrin molecule creating an iron-Tf-TfR complex which enters the cell by endocytosis.

抗原:Synthetic peptide of human TFRC

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: