

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ABCC10 Polyclonal Antibody

多克隆抗体

产品货号: PA08380

第三版

描述:The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, and White). This ABC full-transporter is a member of the MRP subfamily which is involved in multi-drug resistance. Multiple transcript variants encoding different isoforms have been found for this gene.

抗原:Synthetic peptide of human ABCC10

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Isoform 1 is specifically expressed in spleen. Isoform 2 is more widely expressed.