

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PTPN20 Polyclonal Antibody

多克隆抗体

产品货号: PA08385

第三版

**描述:**The product of this gene belongs to the family of classical tyrosine-specific protein tyrosine phosphatases. Many protein tyrosine phosphatases have been shown to regulate fundamental cellular processes. The encoded protein appears to be targeted to sites of actin polymerization. A pseudogene of this gene has been defined on chromosome 10. Alternative splicing results in multiple transcript variants.

**抗原:**Synthetic peptide of human PTPN20B

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Present in many cell lines (at protein level). Widely expressed.