

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ABI3BP Polyclonal Antibody

多克隆抗体

产品货号: PA08396

第三版

描述:TARSH gene expression is dramatically induced in mouse embryonic fibroblasts (MEFs) replicative senescence and suppressed in human lung carcinoma specimens and thyroid carcinomas. The experiment observed that ABI3 expression is reduced or lost in most carcinomas but also that there is a positive correlation between ABI3 and ABI3BP expression. And it express high levels in brain, heart, lung, liver, pancreas kidney and placenta.

抗原:Synthetic peptide of human ABI3BP

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in brain, heart, lung, liver, pancreas kidney and placenta.