

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ADAR Polyclonal Antibody

多克隆抗体

产品货号: PA08414

第三版

描述: This gene encodes the enzyme responsible for RNA editing by site-specific deamination of adenosines. This enzyme destabilizes double-stranded RNA through conversion of adenosine to inosine. Mutations in this gene have been associated with dyschromatosis symmetrica hereditaria. Alternative splicing results in multiple transcript variants.

抗原: Synthetic peptide of human ADAR

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitously expressed, highest levels were found in brain and lung.