

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human POMT1 Polyclonal Antibody

多克隆抗体

产品货号: PA08418

第三版

**描述:**The protein encoded by this gene is an O-mannosyltransferase that requires interaction with the product of the POMT2 gene for enzymatic function. The encoded protein is found in the membrane of the endoplasmic reticulum. Defects in this gene are a cause of Walker-Warburg syndrome (WWS) and limb-girdle muscular dystrophy type 2K (LGMD2K). Several transcript variants encoding different isoforms have been found for this gene.

**抗原:**Synthetic peptide of human POMT1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**Widely expressed. Highly expressed in testis, heart and pancreas. Detected at lower levels in kidney, skeletal muscle, brain, placenta, lung and liver.