

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SPIRE2 Polyclonal Antibody

多克隆抗体

产品货号: PA08431

第三版

描述: SPIRE2, is a 714 amino acid protein belonging to the spire family. Acts as an actin nucleation factor, remains associated with the slow-growing pointed end of the new filament. Involved in intracellular vesicle transport along actin fibers, providing a novel link between actin cytoskeleton dynamics and intracellular transport. Required for asymmetric spindle positioning and asymmetric cell division during meiosis. Required for normal formation of the cleavage furrow and for polar body extrusion during female germ cell meiosis.

抗原: Recombinant protein of human SPIRE2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: