

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ALPPL2 Polyclonal Antibody

多克隆抗体

产品货号: PA08457

第三版

描述: There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The product of this gene is a membrane bound glycosylated enzyme, localized to testis, thymus and certain germ cell tumors, that is closely related to both the placental and intestinal forms of alkaline phosphatase.

抗原: Synthetic peptide of human ALPPL2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Trace amounts in the testis and thymus, and in elevated amounts in germ cell tumors.