

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ALPI Polyclonal Antibody

多克隆抗体

产品货号: PA08494

第三版

描述: There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The intestinal alkaline phosphatase gene encodes a digestive brush-border enzyme. This enzyme is upregulated during small intestinal epithelial cell differentiation.

抗原: Recombinant protein of human ALPI

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: