

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SLC45A3 Polyclonal Antibody

多克隆抗体

产品货号: PA08533

第三版

描述:Solute carrier family 45 member 3 (SLC45A3), also known as prostate cancer-associated protein 6 or prostein, is a protein that in humans is encoded by the SLC45A3 gene. SLC45A3 is expressed in a prostate-specific manner by normal tissues and at a significantly lower level in prostate tumor cell lines. Treatment prostate cancer cell lines with androgens upregulates the expression of SLC45A3.

抗原:Synthetic peptide of human SLC45A3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Prostate specific. Expressed in all prostatic glandular cells. Expressed both in normal and cancerous prostates.