

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human P2RX2 Polyclonal Antibody

多克隆抗体

产品货号: PA08589

第三版

描述:The product of this gene belongs to the family of purinoceptors for ATP. This receptor functions as a ligand-gated ion channel. Binding to ATP mediates synaptic transmission between neurons and from neurons to smooth muscle. Multiple transcript variants encoding distinct isoforms have been identified for this gene.

抗原:Synthetic peptide of human P2RX2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: