

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BFAR Polyclonal Antibody

多克隆抗体

产品货号: PA08600

第三版

描述: A novel protein, BAR, for bifunctional apoptosis regulator, contains domains that enable it to interact with components of both major apoptosis pathways, where it negatively regulates apoptotic signaling. Like the other anti-apoptosis proteins Bap31 and FLIP, BAR contains a DED-like domain that is capable of suppressing apoptosis mediated at the receptor level. In addition, BAR contains a domain that also enables it to interact with the mitochondrial Bcl-2 family of proteins. The presence of these various RING, SAM, DED and TM domains suggests that BAR may serve as a scaffold protein that integrates signaling components of the cells apoptosis-regulatory machinery.

抗原: Synthetic peptide of human BFAR

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed highly in brain, moderately in small intestine, weakly in testes and only faintly in liver and skeletal muscle. Not expressed in heart, kidney, lung and spleen.