

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human C12orf75 Polyclonal Antibody

多克隆抗体

产品货号: PA08607

第三版

描述: C12orf75 (chromosome 12 open reading frame 75), also known as AGD3 or OCC1, is a protein that in humans is encoded by the C12orf75 gene. C12orf75 is high expression in placenta, skeletal muscle, kidney and pancreas tissues.

抗原: Synthetic peptide of human C12orf75

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: