

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SH2D3A Polyclonal Antibody

多克隆抗体

产品货号: PA08616

第三版

描述:SH2D3A (SH2 domain containing 3A), also known as novel SH2-containing protein 1, is a 576 amino acid protein that is thought to play a role in JNK activation. SH2D3A interacts with p130 Cas and is found at low levels in fetal kidney, fetal lung, placenta, adult pancreas, kidney and lung. Subject to post-translational phosphorylation on multiple tyrosine residues, SH2D3A contains one Src homology 2 (SH2) domain. SH2 domains bind to tyrosine-phosphorylated regions of target proteins, frequently linking activated growth factors to putative signal transduction proteins.

抗原:Synthetic peptide of human SH2D3A

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Weakly expressed in placenta, fetal kidney, fetal lung, adult pancreas, adult kidney and adult lung.