

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NEK11 Polyclonal Antibody

多克隆抗体

产品货号: PA08639

第三版

**描述:** This gene encodes a member of the never in mitosis gene A family of kinases. The encoded protein localizes to the nucleoli, and may function with NEK2A in the S-phase checkpoint. The encoded protein appears to play roles in DNA replication and response to genotoxic stress. Alternatively spliced transcript variants have been described.

**抗原:** Recombinant protein of human NEK11

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Poorly expressed in cerebellum, trachea, lung, appendix, and uterus.