

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NLRP9 Polyclonal Antibody

多克隆抗体

产品货号: PA08658

第三版

描述:The protein encoded by this gene belongs to the NALP protein family. Members of the NALP protein family typically contain a NACHT domain, a NACHT-associated domain (NAD), a C-terminal leucine-rich repeat (LRR) region, and an N-terminal pyrin domain (PYD). This protein may play a regulatory role in the innate immune system as similar family members belong to the signal-induced multiprotein complex, the inflammasome, that activates the pro-inflammatory caspases, caspase-1 and caspase-5.

抗原:Synthetic peptide of human NLRP9

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in ileum intestinal epithelial cells. Not detected in peripheral blood mononuclear cells (PubMed:28636595). Expressed in cerebral endothelial cells and, at much lower levels, in brain pericytes (PubMed:28432035).