

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MYBBP1A Polyclonal Antibody

多克隆抗体

产品货号: PA08713

第三版

描述: This gene encodes a nucleolar transcriptional regulator that was first identified by its ability to bind specifically to the Myb proto-oncogene protein. The encoded protein is thought to play a role in many cellular processes including response to nucleolar stress, tumor suppression and synthesis of ribosomal DNA. Alternate splicing results in multiple transcript variants.

抗原: Synthetic peptide of human MYBBP1A

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: