

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MRGPRX1 Polyclonal Antibody

多克隆抗体

产品货号: PA08722

第三版

描述:MRGX1 is an opioid receptor that has been reported exclusively in dorsal root ganglion. MRGX1 is probably involved in the function of nociceptive neurons and may regulate nociceptor function and/or development, including the sensation or modulation of pain. It is potently activated by enkephalins including BAM22 (bovine adrenal medulla peptide 22) and BAM (8-22). BAM22 is the most potent compound and evoked a large and dose-dependent release of intracellular calcium in stably transfected cells.

抗原:Synthetic peptide of human MRGPRX1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Uniquely localized in a subset of small dorsal root and trigeminal sensory neurons.1 Publication