

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NEUROG2 Polyclonal Antibody

多克隆抗体

产品货号: PA08725

第三版

描述: This gene encodes a neural-specific basic helix-loop-helix (bHLH) transcription factor that can specify a neuronal fate on ectodermal cells and is expressed in neural progenitor cells within the developing central and peripheral nervous systems. The protein product of this gene also plays a role in the differentiation and survival of midbrain dopaminergic neurons.

抗原: Synthetic peptide of human NEUROG2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: