

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HMGN5 Polyclonal Antibody

多克隆抗体

产品货号: PA08757

第三版

描述: This gene encodes a nuclear protein with similarities to the high mobility group proteins, HMG14 and HMG17, which suggests that this protein may function as a nucleosomal binding and transcriptional activating protein.

抗原: Synthetic peptide of human HMGN5

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitously expressed.