

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ROS1 Polyclonal Antibody

多克隆抗体

产品货号: PA08851

第三版

**描述:** This proto-oncogene, highly-expressed in a variety of tumor cell lines, belongs to the sevenless subfamily of tyrosine kinase insulin receptor genes. The protein encoded by this gene is a type I integral membrane protein with tyrosine kinase activity. The protein may function as a growth or differentiation factor receptor.

**抗原:** Synthetic peptide of human ROS1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expressed in brain. Expression is increased in primary gliomas.