

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NCAPD3 Polyclonal Antibody

多克隆抗体

产品货号: PA08944

第三版

描述:Condensin complexes I and II play essential roles in mitotic chromosome assembly and segregation. Both condensins contain 2 invariant structural maintenance of chromosome (SMC) subunits, SMC2 and SMC4, but they contain different sets of non-SMC subunits. NCAPD3 is 1 of 3 non-SMC subunits that define condensin II.

抗原:Recombinant protein of human NCAPD3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: