

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ITPR3 Polyclonal Antibody

多克隆抗体

产品货号: PA08967

第三版

描述: This gene encodes a receptor for inositol 1,4,5-trisphosphate, a second messenger that mediates the release of intracellular calcium. The receptor contains a calcium channel at the C-terminus and the ligand-binding site at the N-terminus. Knockout studies in mice suggest that type 2 and type 3 inositol 1,4,5-trisphosphate receptors play a key role in exocrine secretion underlying energy metabolism and growth.

抗原: Synthetic peptide of human ITPR3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in intestinal crypt and villus epithelial cells.