



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AADACL4 Polyclonal Antibody

多克隆抗体

产品货号: PA09067

第三版

描述: Arylacetamide deacetylase-like 4 is a protein in humans that is encoded by the AADACL4 gene. It belongs to the 'GDxG' lipolytic enzyme family. The subcellular location of this protein in the membrane, it is single-pass type II membrane protein. It has signal-anchor, transmembrane and transmembrane helix domain. Function as Hydrolase.

抗原: Synthetic peptide of human AADACL4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: